

Supply Healthier Life!





Latmedical, LLC is an industry leader in the Marketing, Distribution and Supply of high-quality medical products.

We are in the business of well-being, health improvement, and disease prevention by providing quality products and supplies to our patients through our umbrella brand dr. goed.

Dedicated to providing a strong product mix, dr.goed provides a complete and innovative product portfolio solution to the growing healthcare needs within the fields of non-invasive orthopedics, sports medicine, bandaging, wound care and compression therapy.

Our business model has proven to be successful in our targeted regions, not only because we care, but also because we have developed dr.goed.

At latmedical we are very proud of what we do, and thrive to deliver only the best to our patients and customers.

Your Well-being is our Commitment to Excellence.



877 LAT MED0 (528.6330) www.drgoed.com www.latmedical.com









PULL-UPS

Pull-ups are permeable, comfortable items that seem and feel like actual underwear. Ideal for personal use on oneself, they can slide off and on, and provide a relaxed and comfy feeling. For day or overnight use, they come in all shapes and sizes, no matter your requirements; pull-ups are for everyone. Different option are: overnight pull-ups, gender-specific pull-ups, and youth pull-ups.







SPECIFICATIONS

Size	Ü	Ü	S.A.P Weight	Absorbency Weight	Package
	(mm)	(g)	(g)	(ml)	Pcs/bag/ctn
M	720*600	56±5	10	900	10pcs/bag18bags/ctn
L	800*700	63±5	12	1100	10pcs/bag18bags/ctn
XL	880*840	75±5 14	14	1300	10pcs/bag18bags/ctn

MATERIALS

Topsheet: Thermal boned hydrophilic N.W

Absorbent Core: fluff pulp mixed with sap wrapped by

tissue paper+blue

ADL Fluff pulp: Weyerhaeuser Fluff Pulp .

BLD, SAN-DIA, SUMITOMO, BASF etc. as request

Blue ADL: Blue core

Backsheet: Printing colorful pattern PE film backsheet with Wetness Indicator Frontal&side tapes: frontal tape with refastenable adhesive tape

Leakguard: SMS Hydrophobic Non-woven anti-leakage

FUNCTION

- Non woven fabric: keep body surface dry and natural.
- Super absorbent core: fluff pulp + SAP.
- ◆ PE film backsheet: particular surface performance to ensure smooth and gentle touch feeling.
- PP frontal tape: adjustable for fastening and re-fastening with lovely patterns.
- **Wetness indictor:** Pattern fades out when diaper is wet tell you the exact time to change.
- 3D leak guard: protect from leaking.



ADULT DIAPERS

Adult diapers are other very absorbent and convenient options for those facing incontinence. This item has a tendency to work better for people taking care of someone else. Their simple reusable strips make it easier to place them on a loved one with the correct comfort. We ensure options that work for almost everybody. For more safety and extra improvement of permeability, a supporter pad may be positioned inside a diaper with tabs





SPECIFICATIONS

Size	e Length*Width	Total Weight	S.A.P Weight	Absorbency Weight	Package
	(mm)	(g)	(g)	(ml)	Pcs/bag/ctn
М	800*650	83±4	7.5	900	10pcs/bag10bags/ctn
L	960*750	93±4	8	1100	10pcs/bag10bags/ctn
XL	990*800	100±4	9	1300	10pcs/bag10bags/ctn

MATERIALS

Topsheet: Thermal boned hydrophilic N.W

Absorbent Core: fluff pulp mixed with sap wrapped by

tissue paper+blue ADL.

Fluff pulp: Weyerhaeuser Fluff Pulp.

BLD, SAN-DIA, SUMITOMO, BASF etc. as request

Blue ADL: Blue core

Backsheet: Printing colorful pattern PE film backsheet with Wetness Indicator Frontal&side tapes: frontal tape

with refastenable adhesive tape

Leakguard: SMS Hydrophobic Non-woven anti-leakage

FUNCTION

- Non woven fabric: keep body surface dry and
- Super absorbent core: fluff pulp + SAP.
- ♣ Blue ADL: blue Acquisition Distribution Layer can absorb and distribute liquid quickly, and effect prevent from rewet.
- PE film backsheet: particular surface performance to ensure smooth and gentle touch feeling.
- PP frontal tape: adjustable for fastening and re-fastening with lovely patterns.
- Wetness indictor: Pattern fades out when diaper is wet tell you the exact time to change.
- 3D leak guard: protect from leaking.





PADS

Pads are another alternative offering sponginess, wellbeing and convenience. Naturally stack into the user personal underwear with an adhesive backup, they are impeccable for the ones who need protection but doesn't like to wear any of the previous options: pull-up or tab-style diaper. For still more protection, there are the belted pads and the pant liners. We provide also heavy-duty pads for bowel incontinency.







SPECIFICATIONS

Product Name	Non Woven Sanitary Napkin		
Туре	Ultra Thin		
Surface	Soft Non Woven		
Time Used	Day/Night/Liner		
Quality Classification	Low quality/Medium quality		
Length	155-420MM (Specific Length Can Be Customized)		
Absorbency	20ml/30ml/10ml/30ml/10ml/40ml/60ml/10ml		
Inner Packing	PE Individually wrapped		
Special Technology	Pure cotton topsheet & Green silver ion chip		

MATERIALS

Non woven / PE perforated film, SAP, Fluff Pulp, Toilet / Airlaid paper, PE film, Release paper etc.

FUNCTION

- Ultra thin design makes ladies more comfortable and free.
- **♦ Locking-water macromolecule film** top-shgranule to prevent back leakage.
- PE perforated can keeps the surface dry.
- Fluff pulp with SAP can absorb the liquid fast and completely.
- New breathable back-sheet effectively let out wet air and keeps your skin cool.
- Packing can be designed with customer's brand name or logo





UNDERPADS

Underpads are an ideal solution against all types of wetness. These rectangular large pads can be adapted to almost any surface, such as beds, wheelchairs, sofas, chairs for protection against incontinence, and they can be used in hospitals and Dental Clinics as well. They can be used even on pets!









SPECIFICATIONS

Color	blue,green,white,customized
Material	PE film+PP non woven fabric+paper+wooden paper,or +SAP
Weight	20-110g/pc
Package	plastic bag + cartoon bag / + plastic woven bag

CAUTION

- This product is for one-time use and cannot be reused.
- This product must not be used after it is moldy, yellow, or expired, and must not be used if the outer packaging is damaged or contaminated.
- 1t should be placed in a dry place, in a cool, fire-proof and moisture-proof manner.





1.877.LAT MEDO (528.6330) info@drgoed.com www.drgoed.com www.latmedical.com 8461 Lake Worth Road, Suite 434, Lake Worth, FI 33467 USA