

SIZE CHART ACCORDING TO SHOE SIZE FOR COMPRESSION HOSIERY



Please read these instructions carefully to determine the correct size for different product styles.
Measurements should always be taken in the following order, and as soon as possible after waking when swelling is at its lowest.

SIZE	MAN	WOMAN
SMALL	34 - 36	32 - 34
MEDIUM	37 - 39	35 - 37
LARGE	40 - 43	38 - 40
X-LARGE	44 - 47	41 - 42

Latmedical, LLC.

8461 Lake Worth Rd, Suite 434,

Lake Worth, FL 33467 USA

telephone: 1-561-962-4343

toll free: 1-877-LAT-MED0 (528-6330)

SIZE CHART ACCORDING TO CIRCUMFERENTIAL MEASUREMENTS FOR BELOW KNEE, THIGH HIGH AND PANTYHOSE SOCKS



Please read these instructions carefully to determine the correct size for different product styles.

Measurements should always be taken in the following order, and as soon as possible after waking when swelling is at its lowest.

	SMALL	MEDIUM	LARGE	X-LARGE	XX-LARGE
ANKLE	6.5" - 8" 16.5 - 20.5cm	8" - 9.5" 20.5 - 24.5cm	9.5" - 11" 24.5 - 28 cm	11" - 12.5" 28 - 32 cm	12.5" - 14" 32 - 35.5 cm
CALF	10.5" - 14.5" 26.5 - 37 cm	11.5" - 15.5" 29 - 39.5 cm	12.5" - 17" 32 - 43 cm	13.5" - 17.5" 34.5 - 44.5 cm	14.5" - 19.5" 37 - 49.5 cm
THIGH	15.5" - 22" 38 - 56 cm	17.5" - 24" 44.5 - 61 cm	19.5" - 26" 49.5 - 66 cm	22" - 28" 56 - 71 cm	26" - 32" 66 - 81 cm
HIPS	Hasta 40" Hasta 102 cm	Hasta 44" Hasta 113 cm	Hasta 48" Hasta 122cm	Hasta 52" Hasta 132 cm	Hasta 56" Hasta 142cm

Latmedical, LLC.

8461 Lake Worth Rd, Suite 434,

Lake Worth, FL 33467 USA

telephone: 1-561-962-4343

toll free: 1-877-LAT-MED0 (528-6330)

SIZE CHART ACCORDING TO CALF AND KNEE MEASUREMENTS FOR BELOW KNEE SOCKS



Please read these instructions carefully to determine the correct size for different product styles.
Measurements should always be taken in the following order, and as soon as possible after waking when swelling is at its lowest.

		Calf Measurement										
Calk lenght (from the floor to bend in knee)	11"	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"
	(27.5cm)	(30cm)	(32.5cm)	(35cm)	(37.5cm)	(40cm)	(42.5cm)	(45cm)	(47.5cm)	(50cm)	(52.5cm)	(55cm)
14" 35.0 cm	S	S	S	S	M	M	L	L	L	XL	XL	XL
15" 37.5 cm	S	S	S	M	M	L	L	L	XL	XL	XXL	XXL
16" 40.0 cm	S	S	M	M	M	L	L	L	XL	XXL	XXL	XXL
17" 42.5 cm	S	M	M	M	M	L	L	XL	XL	XXL	XXL	XXL
18" 45.0 cm	M	M	M	M	L	L	L	XL	XL	XXL	XXL	XXL
19" 47.5 cm	M	M	M	M	L	L	XL	XL	XL	XXL	XXL	XXL

Latmedical, LLC.

8461 Lake Worth Rd, Suite 434,

Lake Worth, FL 33467 USA

telephone: 1-561-962-4343

toll free: 1-877-LAT-MED0 (528-6330)

SIZE CHART FOR ANTIEMBOLISM STOCKINGS THIGH & THIGH WITH BELT



Please read these instructions carefully to determine the correct size for different product styles.

Measurements should always be taken in the following order, and as soon as possible after waking when swelling is at its lowest.

Ankle Circumference	Calf Circumference	Thigh Circumference	Lenght	Size	Size Color Code (Toes)	Lenght Color Code (Superior)
7" - 7 3/4" (17,8cm - 19,7cm)	Hasta 12 1/2" (31,8cm)	17 1/2" - 25" (44,5cm - 63,5cm)	Up to 29" (73,7cm)	Small - Short	Yellow	Yellow
			29" - 31" (73,7 - 78,7cm)	Small - Regular	Yellow	Natural
			above of 31" (78,7cm)	Small - Long	Yellow	Blue
8" - 8 3/4" (20,3cm - 22,2cm)	12 1/2" - 14 1/2" (31,8cm - 36,8cm)	21" - 28" (53,3cm - 71,1cm)	Up to 29" (73,7cm)	Medium - Short	Natural	Yellow
			29" - 31" (73,7 - 78,7cm)	Medium - Regular	Natural	Natural
			above of 31" (78,7cm)	Medium - Long	Natural	Azul
9" - 9 3/4" (22,9cm - 24,8cm)	14 1/2" - 16 1/2" (36,8cm - 41,9cm)	24" - 31 1/2" (61cm - 80 cm)	Up to 29" (73,7cm)	Large - Short	Blue	Yellow
			29" - 31" (73,7 - 78,7cm)	Large - Regular	Blue	Natural
			above of 31" (78,7cm)	Large - Long	Blue	Blue
10" - 10 3/4" (25,4cm - 27,3cm)	16 1/2" - 18 1/22 (41,9cm - 47cm)	27 1/2" - 33" (69,9cm - 83,8cm)	Up to 31" (78,7cm)	X-large - Regular	Green	Natural
			above of 31" (78,7cm)	X-large - Long	Green	Blue

Latmedical, LLC.

8461 Lake Worth Rd, Suite 434,

Lake Worth, FL 33467 USA

telephone: 1-561-962-4343

toll free: 1-877-LAT-MED0 (528-6330)